


JOYTECH

# Infrared Thermometer

Safe and Accurate



## DET-306( non-contact )



Backlight



1 second reading



Buzzer



10 reading memory



Automatic off

### Product Features

- Measure on forehead
- No probe covers needed
- 10 reading memory
- 1 second reading
- Dual scale
- Compact ergonomic size
- Beeps when ready
- Jumbo display
- Replaceable battery
- Auto-off

### Technical Parameter

Fever Alarm: over 37.8°C(100.4°F)

Unit Weight: Approx. 84.0 grams

Carton Dimension: 32.8x 35.7x 36.5 cm

Range : 34.0°C- 43.0°C(89.6 °F-109.4 °F)

Unit Dimension : 15.4X4.1X4.5cm

Battery : DC3V(2×AAA battery)

Carton Weight: 3.6kgs

Object mode: 0°C- 100°C(32 °F-212 °F)



To measure forehead and object temperature.

